



# Gigogne baby cot

Ref. 1030



Specially designed for use in maternity ward, Médipréma's Gigogne cot can accommodate infants up to the age of 3 months. Médipréma Gigogne baby cot is easy to handle and provides an excellent visibility and ideal access to the baby. Stackable, it saves storage space.

- Detachable transparent PMMA shell with perforated bottom, easy to clean
- Detachable high-resilience comfort mattress : LxWxH=720x340x50 mm
- Bedding height: 767mm
- Continuous head-up/feet-up tilt (+14°/-14°), quickly and gently adjustable using a handle
- White tubular steel frame
- 4 pivoting castors diam.75mm with brakes
- Overall dimensions: LxWxH = 830x551X958mm
- Weigh of entire unit: 14kg
- Maximum permitted load:20kg
- Class I medical device



\*