



Prestige Baby cot

Ref. 1016



Specially designed for use in maternity ward, médipréma's Prestige baby cot can accommodate infants up to the age of 3 months. Médipréma Prestige baby cot is easy to handle and provides an excellent visibility and ideal access to the baby.

- Removable shell made of transparent altuglas, with perforated bottom, easy to clean
- Detachable high-resilience comfort mattress LxWxH=720x340x50 mm
- Bedding height: 686 mm
- Head-up/feet-up tilt (+15°/-18°), quickly and gently adjustable using a handle
- Tubular steel frame
- 4 pivoting castors diam.50mm with brakes
- Overall dimensions: LxWxH = 830x458x877mm
- Weigh of entire unit: 12kg
- Maximum permitted load:20kg
- Class I medical device

