

# Phototherapy masks EyeMax

Ref. 4165 - 4166- 4167



## Comfort and holding quality

The EyeMax phototherapy masks have been designed in order to provide the optimal protection of baby's eyes when it is under phototherapy treatment. The single-piece, wrap-around design attaches at two points on either side of the baby's head and can be independently adjusted to prevent unwanted movement and allow a perfect fit. The special UV light eye protection has been ergonomically designed to prevent any "light leaks".

Available in 3 different sizes in order to adapt to the baby's changing head shape.



	Head measurement
EyeMax <b>Micro</b> ref. 4165	20-25cm
EyeMax <b>Premie</b> ref. 4166	26-32cm
EyeMax <b>Regular</b> ref. 4167	33-38cm

Pack of 20 pieces

Single patient use only – Latex-Free

Components: band (2.449 osy SBL), Eye Pad (Polypropylene Double Sided Cotton Knitted Fabric), Velcro (Nylon hook/Polypropylene Liner)

Maxtec products, distributed by: