

Air/Oxygen Blender MAXBLEND 2

with built-in FiO2 analyzer

Ref : 5111



The MAXBLEND2 is a compact air/oxygen blender integrating a FiO2 analyzer with high and low alarms. It delivers an accurate air/oxygen blend with a large range of flows from 0 to 15L/min. A LCD display directly shows the level of oxygen concentration.

DESCRIPTION

- > **Provided with an acrylic flowmeter with dual scale**
 - 0.25 LPM increments up to 3.0 LPM
 - 1.0 LPM increments to 5 to 15 LPM
- > **Includes Air & O₂ hoses**
- > **Good visibility of the parameters on the display**
 - *Display range* : 0-100% oxygen
 - *Display resolution* : 0.1% oxygen
- > **High and low alarms** : 15-99%
- > **Low battery indicator**

SPECIFICATIONS

Flow range (L/min)	0 to 15
Response time	< 15 seconds at 90% of the final value (at 25°C)
Accuracy	± 3%
Battery life	5000 hours (for continuous use, without alarms)
Power Supply	4x 1.5V batteries
Operating temperature	15°C - 40°C
Weight	Approx. 2.5 kg
Certification	Class IIb medical device

A Maxtec product distributed by :