



For baby with medical gas supplying

Placed in an incubator or an infant warmer, the oxygen hood allows supplying newborns with air and oxygen.

The oxygen hood is made of a transparent clear Plexiglas and grants a good visibility of the baby. Its rectangular shape has smooth edges to prevent baby injury.

4 protected holes allow passing probes and hoses for baby with gas supplying.

Mediprema oxygen hood is available in 2 sizes:

Specifications :

Patient	Size in cm H x L x W	Reference
Premature baby	Small 13 x 17 x 18 cm (3610 cm ³)	3105
Newborn	Medium 15 x 19 x 19.5 cm (5305 cm ³)	3106