



SPECIFICATIONS

médipréma

BIRTHING BED **NIXIE**
REF. 2100

NIXIE, Médipréma's birthing bed, is the ideal woman labor support tool.

Its original "chair and stool" structure offers a unique level of adaptability, perfectly suited to the different positions used during labor phases right up to fetal expulsion (gynaecological position, sideways, on all fours, squatting suspended, etc...).

The **NIXIE** bed enables the woman to be comfortably positioned and gives her a huge freedom of movement, without compromising her comfort and safety. Spacious yet compact, the bed is highly mobile in order to meet any critical emergency situation.

Médipréma proposes in option, a Labor hanger bar, which can be installed without having to rearrange the room. It is detachable and fits directly to the bed structure. It provides a genuine suspension on or near the bed.



BIRTHING BED **NIXIE**

REF. **2100**

DESCRIPTION

The Nixie birthing bed features a chair and a stool.

Electrically actuated:

- Variable chair height : 645-1045mm
- Variable stool height : 500-900mm
- Adjustable height between chair and stool : 0-545mm
- Adjustable backrest: $-12/+70^{\circ}$ (Quick manual Trendelenburg emergency position with release red levers at bed head).
- Adjustable seating: $0/+20^{\circ}$
- Manual tilt on stool: $\pm 15^{\circ}$

The **foot rest stool** is released by a foot-operated pedal, and is just as easy to reconnect. When connected to the chair, it is operated through the general bed remote control. When disconnected, it has its own electrical height adjustment control (optional built-in battery – supplied with wall charger).

The **stirrups** are inserted into an ergonomic and orientable holder attached to the bed. Adjustable height and orientation.

10 liters detachable sliding ABS **placenta tray**.

Retractable, detachable **side safety barriers**, fitted with handles and support zones.

Standardized side rails on the seat and stool, to which numerous accessories can be attached.

Wired remote control. Detachable, rail mounted flexible support, for easy attachment and directional positioning.

125mm diameter double roller castors with **central braking** by side pedal on the chair. 100mm diameter castors with brake on the stool.



Pair of grab handles with quick-release multi-directional clamps, to set on the side rails of either the seat or the stool. Uniquely versatile to use, due to their shape and endless adjustment positions: grab handles, foot/heel rests, or used as handles to move the bed easily around.

Comfortable, high-ductility polyurethane **foam mattress**, 8cm thick. 2 separate sections for chair and stool. Antibacterial polyurethane cover with welded seams. Fireproof. Small indentation for the pelvis area.

Battery backup with low battery alarm.

Overall dimensions with barriers:

(WxL) 2080mm x 1025mm

Bed size: 2040mm x 860mm

Chair dimensions: (WxL) 1525mm x 1025mm

Maximum load: 200kg

Electricity supply: 230V-50/60Hz (optional 110V)

Electrical classification: class 1

Class 1 medical device

OPTIONS

- Independent stool with electrical height adjustment (built-in battery - supplied with wall charger)
- Stainless steel telescopic IV pole - 4 hooks - and base attachment. Can be freely positioned and locked on either side of the bed head.
- Footrest for forceps
- 2nd remote control

ACCESSORIES

- Detachable labour hanger bar. Supplied with suspension sling and detachable central strap.
- Handrail (elbow high) on quick-release multi-directional clamps
- Footrests
- Additional pair of grab handles on quick-release multi-directional clamps
- Flexible support for the wired remote control
- Strap for stirrups
- Wall mount support for accessories
- Positioning cushion



BIRTHING BED **NIXIE**

REF. **2100**

POSITIONS OF LABOR AND DELIVERY WITH NIXIE

According to the APNOR® method of De Gasquet (Postural-Respiratory Method)

Mobility



On all fours and squatting



Suspended postures



Sideways



“Adapted” Gynaecological position

