médipréma

Specifications

FOAM & MICRO-**GRANULES TRANSPORT MATTRESS**

Réf. 5044




Support and comfort for newborn babies

By combining both a comfortable shape recall foam mattress with polystyrene microgranules acknowledged for their maintaining properties, the mattress Carewave by Médipréma brings an effective solution to the problems of safety and comfort during transport.

- The shape recall mattress is made of visco-elastic foam which moulds the shape of the child's body. Its soft comfort enables to reduce inconvenience caused by jolts during transport.
- The additional micro-granules filling enables to maintain the newborn body. Better retained the baby feels more secure and more quiet during transport.

Specially sized to be used in médipréma Nite transport incubator, the foam & micro-granules transport mattress. placed directly on the mattress, provides ideal comfort conditions to receive the baby.



médipréma



The foam and micro-granules transport mattress can be attached to the bedding area thanks to the Magic Strips straps. The Magic Strips is a limp micro-fiber velcro strap, soft to the touch, very easy to put in place as well as to take off.

The foam and micro-granules transport mattress is **light** and **easy to put in place**. It can be taken off as easily in order to quickly install the baby on the hard bedding for special medical procedures.

The bi-elastic polyurethane cover is **suited to skin contact**, and is **easy to clean** and disinfect by sponging.

Product	Specifications	Reference	Packing
Foam & microgranules transport mattress CareWave médipréma	 Dimensions: 61 x 38 cm Visco-elastic polyurethane foam mattress Polystyrene micro-granules Bi-elastic polyurethane cover Water resistant Double seam Cleaning by sponging 	• 5044	• Per unit

