



# Thermal reflector probe patch

Ref : 3159

### Single use probe patch

#### Adhesive hydrogel

- > Secures temperature probe to the skin
- > Does not disrupt skin integrity when removed
- > Will not leave residue on the skin
- > Able to be repositioned

#### Thermal reflector

- > Gold reflective foil reflects heat for accurate temperature monitoring



### Description :

Composition - Above	Reflective polyester film
Composition - Below	Adhesive hydrogel
Latex free	Yes

Probe patch Diameter	3 cm
Packaging	Individually pouched Packed by 200 units
Certification	CE Marked



An Axelgaard Manufacturing Co,Ltd. product distributed by :



ZA Node Park Touraine - 470 rue Gilles de Genes - BP 8 - 37310 TAUXIGNY - FRANCE  
 Tél +33 (0)2 47 28 47 87 - Fax +33 (0)2 47 27 35 85 - www.mediprema.com