

# Calibed

Warming mattress



Calibed warming mattress provides the ideal and constant temperature that newborn babies need to feel secure. This protective heat also offers feather-soft comfort for baby. The Calibed warming mattress is removable and can be used in any type of cradle. Calibed provides child's skin temperature monitoring. The high and low alarms prevent any risk of hypothermia or hyperthermia.

# Calibed - Ref.1085

Warming mattress



## SPECIFICATIONS

**Mattress temperature regulation** from 28 to 38.5°C  
(with secure access under 35°C and above 37°C)

**Uniform warming** over the entire bedding area

**Skin temperature monitoring**

**Sound and visual alarms.** Advanced, offset management system for alarms (Quietude mode). 3 sets of visual and sound alarms have been defined depending on the level of emergency for the comfort of newborn baby and the medical team.

**Heating power :** 75 watts (absorbed power : 120 VA)

**Power supply :**

- Control box : 220-240V / 50-60Hz
- Mattress : 24VAC security low voltage

**Electrical Classification :** class I

**Class IIb medical device**



Calibed positionned in a care crib (ref.1020)

## DESCRIPTION

**Removable solution.** The warming mattress is made up of a control box, a warming element and a gel mattress.

**Control box with LCD screen** for a good visibility of data and alarms. A handle provides to fix on any type of cradle.

**Very comfortable polyurethane gel mattress** (latex- and silicone-free)

**Delivered with a removable cocoon** including a cotton cover and a microgranules flexitube cushion. The baby is comfortably installed in a safe environment.

**Dimensions :**

- Control box without handle (L x W x H) : 300 x 190 x 75 mm
- Heating element : 600 x 330 x 10 mm
- Gel mattress : 615 x 385 x 13 mm

**Delivered with one skin probe**



Calibed positionned in a baby cot (ref.1016 or 1030)  
with included cocoon for baby's comfort